

HS-1500 Automatic Tissue Processor



Features:

1. English menu, man-machine interaction, simple operation.
2. Prompt judge and automatic processing in case of abnormal condition, which can prevent the tissue from being damaged.
3. Double-control heat without water
4. Based on different tissues, several running programs available for customers
5. Self-provided power supply can continue to operate for 4 to 6 hours during a power outage. The tissue protection cylinder can be automatically operated after a power outage. Does not contain heated paraffin. If the temperature of the heated wax cylinder does not reach the set temperature, the tissue will not enter the wax cylinder and return to the protection cylinder to wait. High degree of intelligence
6. Environment friendly: The air as purified through internal circulation can be exhausted off the window.
7. With transmission of gear and rack, mechanical failures can be reduced.
8. Several standby programs available

Main Technical Data:

1. Tanks : 12 Tanks (3 wax Tanks in the end)

2. Capacity of single tank : 1500ml
3. Basket can hold 75 specimens
4. Single-tanks processing time: Randomly set within the range from 0 to 99 hour(s).
5. The temperature of wax cylinder can be randomly set within the range of 99°C.
6. Time for cylinder to be lifted: 38 seconds
7. Tissue-protection cylinder sequence: 7th cylinder
8. Times of mixing per minute: 3 times

Accessories:

Standard Setting List

- | | |
|---------------------|--------|
| 1. Tanks | 12 pcs |
| 2. Tissue basket | 1 pc |
| 3. The exhaust pipe | 1 pc |
| 4. Tissue cassette | 20 pcs |
| 5. Fuse (5 A) | 2 pcs |
| 6. Power line | 1 line |

Optional Accessories

1. Extra Tank
2. Tissue Basket

Net Weight: 67kg

Dimensions : 1130*400*500mm

Package Info:

Each instrument packed Wooden Case

Gross Weight: 100 kgs

Dimensions: 1260*570*640 mm,